

Sensacare Product Specification Sheet

SENSAWATCH CONNECT



Product Overview

The SensaWatch Connect is a professional-grade 4G smartwatch designed for healthcare monitoring. It supports fingertip SpO2 tracking, GPS positioning, fall alerts, and various health vitals including heart rate, sleep, and body temperature monitoring, all in a durable waterproof form factor.



Product Name:
SensaWatch Connect
Model: SW-C102
Category: 4G Fingertip
Health Smartwatch



Display: 1.83 inch IPS,
240x280 px IPS, Goodix
3018 heart rate sensor
Battery: 380mAh with 3-7
days usage




Charging: Magnetic
charging, 2 hours

Bluetooth: BLE 5.2

Standby Time: 30 days

Key Features

- Fingertip SpO2 measurement
- Blood pressure, temperature, heart rate, sleep
- Built-in GPS + LBS + WiFi positioning
- SOS call, geofencing, and fall detection
- Two-way voice communication
- Nano SIM support with GSM + 4G LTE
- Compatible with Bluetooth accessories
- Waterproof IP67



Various Bluetooth accessories can be connected via Bluetooth, including a blood pressure monitor, glucose monitor, body weight scale, and pulse oximeter.



Package Contents:
SensaWatch Connect smartwatch, charging cable, user manual, gift box.
LDS+PC+ABS watch case, silicon strap.



Certifications & Compliance:
CE/ROHS/FCC Bluetooth
Certified with a 1-year standard warranty.

Environmental Conditions

The watch features an IP67 waterproof rating for daily splash resistance and wear, and is constructed with an LDS+PC+ABS case and a comfortable silicon strap.