

Sensacare Product Specification Sheet

SENSASCOPE STETHOSCOPE



Product Overview

The SensaScope Stethoscope enables professional heart and lung auscultation at home. It connects wirelessly to mobile devices for real-time listening and data sharing, with built-in ambient noise reduction and auto recording features.



Product Name:
SensaScope Stethoscope
Model: SBS-21 Category:
Digital Wireless
Stethoscope



Size: 45 x 45 x 26 mm,
Weight: 78g, Battery:
200mAh, Interface:
Bluetooth and 3.5mm
headphone/recording
ports 4-7 days usage
Frequency Range: 20Hz –
3kHz




Charging: Magnetic USB,
5V/500mA
charging, 2 hours
Bluetooth: BLE 5.0

Standby Time: 40 days

Key Features

- Wireless Bluetooth auto connectivity
- Automatic recording when contact is detected
- Noise detection alerts for cleaner recordings
- Telemedicine-ready with remote data sharing
- One-button operation (power/mode switch)
- Ambient Noise Reduction (ANR)
- Up to 24x sound amplification
- Dual-mode auscultation: Diaphragm / Bell
- Suitable for: Adult, Infant, Pediatric
- Noise Elimination: Up to 85%



Chestpiece: Kirsite,
Polyurethane-coated
silicone diaphragm
Frequency Range: 20Hz –
3kHz



Package Contents:
SensaScope Stethoscope,
charging cable, user
manual, Optional 3.5mm
Headphones and Recording
Cable



Certifications & Compliance:
CE/ROHS/FCC Bluetooth
Certified with a 2-year
standard warranty.

Environmental Conditions

Operating Temp: -10°C to 40°C (14-104 °F) | Humidity: 5-95% RH, non-condensing
Storage Temp: -25°C to 70°C