

Sensacare Product Specification Sheet

SENSAWATCH BALANCE



Product Overview

The SensaWatch Balance is a modern health monitoring smartwatch designed for fitness enthusiasts and wellness-conscious users. It offers a variety of biometric tracking capabilities in a compact, stylish form.



Product Name:
SensaWatch Balance
Model: SW-B600

Category: Health
Monitoring Smart Watch



1.47" HD display, 172 x
320px, full touch screen,
160mAh battery

30 day standby, 4-7 days




Magnetic pin charging
Charge time: 2 hours
Bluetooth: BLE 5.0

Standby Time: 40 days

Key Features

- Measures blood glucose, ECG, HRV
- Body temperature, blood oxygen, HR
- Blood pressure, sleep, female health
- 15 sports modes
- Smart notifications
- Call reminder, message push, remote camera



Supports multi-languages



Package Contents:
SensaWatch Balance, User
Manual, Charging Cable +
Gift Box



Certifications & Compliance:
CE, ROHS.

1-year standard warranty.

Environmental Conditions

The watch features an IP67 waterproof rating for daily splash resistance and wear, and is constructed with an LDS+PC+ABS case and a comfortable silicon strap.