

Sensacare Product Specification Sheet

SENSAWATCH MOVE



Product Overview

The SensaWatch Move is a stylish and compact health-focused smartwatch designed to support wellness tracking. It includes key biometrics like heart rate and oxygen saturation while offering smartwatch functionality such as notifications and app connectivity.



Product Name:
SensaWatch Move
Model: SW-M420
Category: Blood sugar ECG
blood oxygen health smart
watch



Display: 1.39" HD Display,
360x360 px, Full Touch
Screen, Battery: 280mAh
with 30 days standby, 4-7
days usage



Charging: Magnetic pin
charging, 2 hours
Bluetooth: BLE 5.0

Standby Time: 40 days

Key Features

- Health data: blood glucose + ECG + HRV
 - Body temperature + blood oxygen + heart rate
 - Blood pressure + sleep + female health
 - Non-invasive blood glucose
 - ECG measurement: wrist, chest, heart rate belt
- Osram SFH2201 + TYHX3220 sensor
 - Call reminder, message push, remote camera
 - Full touch + side button controls
 - Standby up to 40 days
 - 10 sports modes



Multilingual Support: 14
device languages and 19
app languages including
English, Italian, Spanish,
French, German,
Portuguese, Chinese,
Japanese, and more.



Package Contents:
SensaWatch Move
smartwatch, charging cable,
user manual, and gift box.
Materials include alloy case
with TPU/leather strap
options.



Certifications & Compliance:
CE/ROHS/FCC Bluetooth
Certified with a 1-year
standard warranty.

Environmental Conditions

This device meets the IP67 standard, ensuring it is fully protected against dust ingress and can withstand immersion in water up to 1 meter for 30 minutes, ideal for daily wear.