



# Sensacare RPM: Transforming Remote Patient Monitoring

Smarter monitoring. Better care, anywhere.

**Sensacare** is a flexible, real-time health monitoring system designed for care providers, clinics, and health services. It brings together device-agnostic tracking, automated insights, and tailored care planning into a single, easy-to-use platform.

# Key Features



## Real-Time Vitals Tracking

Monitors 13+ vitals: heart rate, SpO<sub>2</sub>, temperature, blood pressure, respiration, HRV, and more

Works with Apple Watch, Samsung, BLE wearables and medical-grade sensors

Automatically adjusts to connected devices



## Health Trends and Insights

Highlights changes before they become serious

Predicts risk across chronic conditions

Identifies potential issues using built-in trend analysis



## Built-In Video Consults

One-click telehealth built into the platform

Share notes, review vitals, and consult remotely

Supports calls with multiple care team members



## Custom Alerts and Safety Thresholds

Set personalised limits for each client

Smart alert logic reduces unnecessary alarms

Immediate notifications via email, SMS, or dashboard



## Dashboard Heatmaps and Trends

Spot vital changes at a glance

View historical patterns and performance

Easy-to-read sparkline charts and visual cues



## Care Plans and Reporting

Create care plans tailored to individual needs

Weekly and monthly summaries

Reports are ready to print or upload to EHRs

# More Features

## User Access Management

- **Carers:** live alerts, assigned clients, latest vitals
- **Admins:** oversight of all data, export tools, full access
- **Super Admins:** organisation-wide settings and analytics

## Offline Functionality

- Records data without internet
- Syncs automatically when connection is restored
- Ideal for mobile workers and rural care

## Security and Compliance

- Compliant with HIPAA and Australian Privacy Principles
- Access is role-restricted and fully tracked
- All data is encrypted

## System Integrations

- Supports FHIR and major EHR systems
- Open API for device makers and health tech partners
- Integrates with CDMP and aged care workflows

## Device Compatibility

- Optimised for professional use
- Works in the browser – no app needed
- Touch-ready for tablets used in the field

# Real-Time Vitals Tracking



## Comprehensive Monitoring

Tracks 13+ vitals including heart rate, SpO<sub>2</sub>, and blood pressure.



## Device Flexibility

Works with Apple Watch, Samsung, and medical-grade sensors.



## Adaptive Technology

Automatically adjusts to connected devices with no manual setup.



# Custom Alerts and Thresholds

Sensacare allows care teams to set **custom threshold ranges for each individual**, rather than relying on one-size-fits-all presets. This means alerts are tailored to the patient's unique health profile, whether it's adjusting heart rate parameters for someone with a cardiac condition or setting temperature thresholds for infection-prone clients.



## Personalised Limits

Set custom thresholds for each patient's unique health profile.




## Smart Alert Logic

Advanced algorithms reduce false alarms and alert fatigue.



## Immediate Notifications

Receive critical alerts via email, SMS, or dashboard.



# Health Trends and Insights



## Early Detection

Highlights subtle changes before they become serious complications.

Provides early intervention opportunities to prevent emergencies.



## Risk Prediction

Uses AI to predict potential risks across chronic conditions.

Identifies patterns that might be missed in routine monitoring.



## Trend Analysis

Built-in analytics identify concerning patterns automatically.

Helps prioritize patients needing immediate attention.

# Built-In Video Consultations



## One-Click Telehealth with Live Vitals

Start secure video consultations directly from the platform with live vitals presented for that client.

## Collaborative Care

Share notes, review vitals, and consult remotely with the care team.

## Multi-Provider Support

Facilitate calls with multiple care team members for complex cases.



# Competitor Snapshot

Feature	Sensacare	Withings RPM	CareMonitor	Telstra Health	Cardiab / ResApp
Real-Time Vitals Monitoring	✓	Limited vitals	Basic vitals only	Basic vitals only	Condition-specific
Custom Alerts & Thresholds	✓	✗	Preset only	Preset only	✗
AI Health Insights	✓	✗	✗	✗	✗
Telehealth Integration	✓	✗	Limited	Optional	✗
Care Plan Automation	✓	✗	✗	✗	✗
Device Agnostic	✓	✗	✗	✗	✗
EHR/API Integration	✓	✗	Limited	Limited	✗
Offline Sync	✓	✗	✗	✗	✗
Training CMS	✓	✗	✗	✗	✗
Chronic Disease Tools	✓	✗	✗	✗	✓ Cardiac focus
Zero Upfront Cost Option	✓	✗	✗	✗	✗
Compliance (HIPAA, APP)	✓	✓ HIPAA only	✓	✓	✓

Sensacare stands well ahead of its competitors by offering a more comprehensive, flexible, and future-ready remote monitoring solution. While many platforms provide only basic vitals or fixed alert systems, Sensacare delivers dynamic, real-time health tracking with personalised thresholds, AI-driven insights, built-in telehealth, and seamless EHR integration, all without locking users into specific devices. With offline sync, zero upfront cost options, and tools to automate chronic disease care, Sensacare goes beyond monitoring to actively support better outcomes for patients and providers alike.



# Role Based

The Sensacare platform is built with a **role-based architecture**, ensuring each user sees only the information relevant to their responsibilities. Carers are provided with a clear dashboard showing assigned clients, real-time alerts, and vital signs that require attention. Admins have broader oversight, including full client management, reporting tools, and access to audit logs. A dedicated **Super Admin role** (coming soon) will offer organisation-level controls, including analytics, system configuration and role permissions across multiple teams or facilities. This structure improves data security, simplifies user experience, and ensures accountability at every level of care delivery.



## Carers

See live alerts, assigned clients, and latest vitals.



## Administrators

Access oversight of all data, export tools, and full system control.



## Super Admins

Manage organization-wide settings and analytics (coming soon).

# Seamless Security and Connectivity



## Secure and Compliant

Meets HIPAA and Australian Privacy Principles standards.



## Works Offline


Records data without internet and syncs when connection returns.



## Open Integrations

Supports FHIR, major EHR systems, and offers open APIs.

# Why Healthcare Teams Choose SensaCare



# Training

Sensacare includes a built-in **training and onboarding module** designed to get care teams up and running quickly. Whether you're a first-time user or managing a large rollout across multiple sites, the platform guides users through key features, workflows and compliance requirements. Interactive tutorials, step-by-step walkthroughs, and role-specific instructions ensure carers, clinicians and admins gain the confidence to use the system effectively. This in-platform training reduces the need for external support, shortens onboarding time, and ensures consistent standards of care across your organisation.



# Contact Us

We're here to help with all your remote patient monitoring needs.

## Visit Our Website

[www.sensacare.com.au](http://www.sensacare.com.au)

## Get a Demo

[Schedule a personalised demonstration](#)

## Join our Pilot Program

[Experience Sensacare firsthand](#)

## Contact Us

[www.sensacare.com.au/contact](http://www.sensacare.com.au/contact)

Available Monday-Friday, 9am-5pm